



DEMULSIFICATION TEST APPARATUS FOR LUBRICATING OILS

APPLICATION

Tests the ability of medium to high viscosity oils to separate from water, in presence of contamination and turbulence.

COMPLIANCE : ASTM D 2711, DIN 51353

FEATURES

Variable speed stirrers with Digital RPM Indicators


All wetted parts made of SS

Full view digital constant temperature bath with sample holders.



STANDARD ACCESSORIES

- 11020-00 Constant temp. bath. 230V operation.
- 11020-02 Glass Separating funnels with provision for aligned suspension.
- 11020-03 High speed sample stirrers with variable speed.



Weight :

Size :

OPTIONAL ACCESSORIES

- 11020-01 Constant temp. bath. 110V operation
- 11020-04 Centering Gauge
- 11020-05 Solvent tank with forced air drier for cleaning the stirrer assembly.
- 11020-06 Constant temp. conditioning bath 2/4/6/8 vessel
- xxx High Speed centrifuge
- xxx Graduated centrifuge tube
- xxx Thermometer ASTM9F
- xxx Thermometer ASTM9C
- xxx Digital tachometer non contact type



Culture Instruments India LLP

37, III MAIN ROAD, VYALIKAVALL, BANGALORE - 560 003. (INDIA)

PH: 91 - 080 - 23349239 / 23449007 / 09886780513

Email: cultureinstruments@gmail.com / sales@technoinstruments.co.in

Website: www.technoinstruments.co.in

Works: D-428, 10th MAIN, PEENYA 2nd STAGE, P.I.A, BANGALORE - 560 058. (INDIA)

* Products & Design Subject To Change Without Prior Notice.